


Anti-Abi1 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q8Izp0

Item Number: CM0021121



Introduction

Affinity-purified rabbit polyclonal antibody targeting human ABI1 (Abi-1, E3B1) for WB, IHC-P, ELISA. Detects endogenous ABI1 across human, mouse, and rat samples—suitable for studies on cell signaling, cytoskeletal dynamics, and neurodevelopment.

[Learn More](#)

Parameter	Value
Product Name	ABI1 Rabbit pAb
Remarks/Alias	E3B1; ABI-1; ABLBP4; NAP1BP; SSH3BP; SSH3BP1; ABI1
Species	Human
Gene ID (Human)	10006
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 180-230 of human ABI1 (NP_001012770)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8IZP0
Protein Weight	55kDa
Shipping	Ice bag