

Anti-Abi3 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q9P2A4

Item Number: CM0006602

Introduction

Rabbit monoclonal antibody targeting human ABI3 (NESH/SSH3BP3), suitable for Western blot, immunohistochemistry (paraffin), and ELISA. Detects human, mouse, and rat ABI3 with predicted molecular weight 39kDa. Ideal for tumor metastasis and cell motility studies.

[Learn More](#)

Parameter	Value
Product Name	ABI3 Rabbit mAb
Remarks/Alias	NESH; SSH3BP3; ABI3
Species	Human
GeneID (human)	51225
GeneID	51225
Immunogen	Recombinant fusion protein containing a.a. 197-366 of human ABI3 (NP_057512.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9P2A4
Protein Weight	39kDa
Shipping	Ice bag