


Anti-Abl2 Rabbit Monoclonal Antibody For Wb, Elisa - P42684

Item Number: CM0002145



Introduction

High-quality rabbit monoclonal antibody targeting human ABL2 (ARG) for WB and ELISA applications. Cross-reacts with rat. Validated for cytoskeleton remodeling and cell adhesion research.

[Learn More](#)

Parameter	Value
Product Name	ABL2 Rabbit mAb
Remarks/Alias	ARG; ABLL; ABL2
Species	Human
GeneID (Human)	27
GeneID	27
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 500-600 of human ABL2 (P42684).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P42684
Protein Weight	128kDa
Shipping	Ice bag