

## Anti-Ablim2 Rabbit Polyclonal Antibody For Wb, Elisa - Q6H8Q1

Item Number: CM0021795



### Introduction

High-quality anti-ABLIM2 rabbit polyclonal antibody validated for Western blot (WB) and ELISA. Reacts with human, mouse, and rat samples. Ideal for skeletal muscle research and actin-binding LIM protein 2 detection.

[Learn More](#)

| Parameter        | Value                 |
|------------------|-----------------------|
| Product Name     | ABLIM2 Rabbit pAb     |
| Remarks/Alias    | ABLIM2                |
| Species          | Human                 |
| GeneID (Human)   | 84448                 |
| GeneID           | 84448                 |
| Immunogen        | Recombinant protein   |
| Source           | Rabbit                |
| Category         | Polyclonal Antibodies |
| Application      | WB, ELISA             |
| Cross Reactivity | Mouse, Rat            |
| SWISS            | Q6H8Q1                |
| Protein Weight   | 68kDa                 |
| Shipping         | Ice bag               |