


Anti-Ablim3 Polyclonal Antibody For Wb And Elisa - 094929

Item Number: CM0020395



Introduction

Rabbit polyclonal antibody targeting human ABLIM3, an actin-binding LIM scaffold protein. Applications: WB, ELISA. Cross-reacts with mouse, rat. Ideal for cardiac and brain research. Offered by CamelBio.

[Learn More](#)

Parameter	Value
Product Name	ABLIM3 Rabbit pAb
Remarks/Alias	HMFN1661; ABLIM3
Species	Human
Gene ID (Human)	22885
Gene ID	22885
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 350-520 of human ABLIM3 (NP_001287957.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	094929
Protein Weight	78kDa
Shipping	Ice bag