

Anti-Abr Rabbit Polyclonal Antibody For Wb And Elisa - Q12979

Item Number: CM0022161

Introduction

High-quality rabbit polyclonal antibody against human ABR (Active breakpoint cluster region-related protein). Validated for WB and ELISA, cross-reacts with mouse and rat. Ideal for GTPase signaling and neuronal research.

[Learn More](#)

Parameter	Value
Product Name	ABR Rabbit pAb
Remarks/Alias	MDB; ABR
Species	Human
GeneID	29
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 200 to the C-terminus of human ABR (NP_001243776.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q12979
Protein Weight	98kDa
Shipping	Ice bag