

Anti-Abtb1 Rabbit Polyclonal Antibody For Wb, Elisa - Q969K4

Item Number: CM0024472

Introduction

High-quality Anti-ABTB1 Rabbit Polyclonal Antibody validated for WB and ELISA, with cross-reactivity to mouse and rat. Useful for studying PTEN growth-suppressive signaling and developmental processes.

[Learn More](#)

Parameter	Value
Product Name	ABTB1 Rabbit pAb
Remarks/Alias	BPOZ; BTB3; BTBD21; EF1ABP; PP2259; ABTB1
Species	Human
GeneID (Human)	80325
GeneID	80325
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-140 of human ABTB1 (NP_742024.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q969K4
Protein Weight	54kDa
Shipping	Ice bag