

Anti-Acaa1 Polyclonal Antibody For Wb, Elisa - P09110

Item Number: CM0022197



Introduction

Rabbit polyclonal antibody targeting human ACAA1 (3-ketoacyl-CoA thiolase), validated for WB and ELISA with cross-reactivity to mouse. Suitable for studying peroxisomal fatty acid beta-oxidation.

[Learn More](#)

Parameter	Value
Product Name	ACAA1 Rabbit pAb
Remarks/Alias	ACAA; THIO; PTHIO; Lnc-Myd88; ACAA1
Species	Human
GeneID (Human)	30
GeneID	30
Immunogen	Synthetic peptide corresponding to amino acids 210-310 of human ACAA1 (NP_001598.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P09110
Protein Weight	44kDa
Shipping	Ice bag