

Anti-Acaa1 Rabbit Polyclonal Antibody For Wb, Elisa - P09110

Item Number: CM0026955



Introduction

ACAA1 Rabbit pAb, validated for WB and ELISA, targets peroxisomal 3-ketoacyl-CoA thiolase in human and rat. Ideal for fatty acid metabolism and peroxisomal beta-oxidation studies.

[Learn More](#)

Parameter	Value
Product Name	ACAA1 Rabbit pAb
Remarks/Alias	ACAA; THIO; PTHIO; Lnc-Myd88; ACAA1
Species	Human
GeneID	30
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-300 of human ACAA1 (NP_001598.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P09110
Protein Weight	44kDa
Shipping	Ice bag