


Anti-Acaa2 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P42765

Item Number: CM0002187



Introduction

High-quality rabbit monoclonal antibody targeting human ACAA2 (3-ketoacyl-CoA thiolase). Validated for WB, IF/ICC, and ELISA; cross-reacts with mouse and rat. Ideal for fatty acid metabolism and mitochondrial research.

[Learn More](#)

Parameter	Value
Product Name	ACAA2 Rabbit mAb
Remarks/Alias	DSAEC; ACAA2
Species	Human
Gene ID (Human)	10449
Gene ID	10449
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 298-397 of human ACAA2 (P42765).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P42765
Protein Weight	42kDa
Shipping	Ice bag