


Anti-Acaa2 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P42765

Item Number: CM0010364



Introduction

Rabbit monoclonal ACAA2 antibody validated for WB, IHC-P, ELISA. Detects human, mouse, rat 3-ketoacyl-CoA thiolase (42kDa), a mitochondrial fatty acid beta-oxidation enzyme. Ideal for IVD and metabolism research.

[Learn More](#)

Parameter	Value
Product Name	ACAA2 Rabbit mAb
Remarks/Alias	DSAEC
Species	Human
Gene ID	10449
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 17-271 of human ACAA2 (NP_006102.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P42765
Protein Weight	42kDa
Shipping	Ice bag