


Anti-Acaa2 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P42765

Item Number: CM0021020



Introduction

High-quality rabbit polyclonal antibody against human ACAA2, validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for mitochondrial beta-oxidation and fatty acid metabolism studies.

[Learn More](#)

Parameter	Value
Product Name	ACAA2 Rabbit pAb
Remarks/Alias	DSAEC; ACAA2
Species	Human
GeneID (Human)	10449
GeneID	10449
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human ACAA2 (NP_006102.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P42765
Protein Weight	42kDa
Shipping	Ice bag