


Anti-Acac1/Acc1 Polyclonal Antibody For Ihc-P And Elisa - Q13085

Item Number: CM0022721



Introduction

Rabbit polyclonal antibody against human ACAC1 (ACC1), validated for IHC-P and ELISA. Cross-reacts with mouse and rat. Suitable for fatty acid biosynthesis and metabolic research applications.

[Learn More](#)

Parameter	Value
Product Name	ACAC1 Rabbit pAb
Remarks/Alias	ACC; ACAC; ACC1; ACCA; Acac1; hACC1; ACACAD; ACCalpha; ACACalpha; acetyl-ACC1
Species	Human
Gene ID (Human)	31
Gene ID	31
Immunogen	Recombinant protein corresponding to a sequence within amino acids 1-75 of human ACAC1 (NP_942133.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13085
Protein Weight	266kDa
Shipping	Ice bag