

Anti-Acad11 Rabbit Polyclonal Antibody For Wb And Elisa - Q709F0

Item Number: CM0013012



Introduction

Rabbit polyclonal antibody targeting human ACAD11, validated for Western blot and ELISA. Cross-reacts with mouse and rat. Detects an ~87 kDa protein involved in CNS fatty acid beta-oxidation. Suitable for metabolic research.

[Learn More](#)

Parameter	Specification
Product Name	ACAD11 Rabbit pAb
Remarks/Alias	ACAD-11; ACAD11
Species	Human
Gene ID (Human)	84129
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-230 of human ACAD11 (NP_115545.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q709F0
Protein Weight	87kDa
Shipping	Ice bag