

Anti-Acad8 Polyclonal Antibody For Wb, Elisa - Q9Z0M6

Item Number: CM0019228

Introduction

High-quality rabbit polyclonal antibody targeting ACAD8 (IBDH/ARC42) for WB and ELISA applications. Exhibits cross-reactivity with mouse and rat. Suitable for research on leukocyte adhesion and migration.

[Learn More](#)

Parameter	Value
Product Name	ACAD8 Rabbit pAb
Remarks/Alias	IBDH; ARC42; ACAD-8; ACAD8
Species	Human
GeneID (human)	27034
GeneID	27034
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-415 of human ACAD8 (NP_055199.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9Z0M6
Protein Weight	45kDa
Shipping	Ice bag