

Anti-Acad8 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q9Uku7

Item Number: CM0004728

Introduction

High-quality rabbit monoclonal antibody against human ACAD8 (isobutyryl-CoA dehydrogenase). Validated for WB, IHC-P, ELISA; cross-reacts with mouse and rat. Ideal for studying valine catabolism and metabolic disease.

[Learn More](#)

Parameter	Value
Product Name	ACAD8 Rabbit mAb
Remarks/Alias	IBDH; ARC42; ACAD-8
Species	Human
GeneID	27034
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-415 of human ACAD8 (NP_055199.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UKU7
Protein Weight	45kDa
Shipping	Ice bag