

Anti-Acad9 Rabbit Polyclonal Antibody For Wb And Elisa - Q9H845

Item Number: CM0026155



Introduction

Rabbit polyclonal antibody targeting human ACAD9, validated for WB and ELISA. Cross-reacts with human and mouse. Essential for mitochondrial complex I assembly and fatty acid oxidation studies.

[Learn More](#)

Parameter	Value
Product Name	ACAD9 Rabbit pAb
Remarks/Alias	NPD002; MC1DN20; ACAD9
Species	Human
Gene ID	28976 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human ACAD9 (NP_054768.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9H845
Protein Weight	69kDa
Shipping	Ice bag