


Anti-Acadl Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P28330

Item Number: CM0016571



Introduction

Rabbit polyclonal antibody against human ACADL (LCAD). Validated for WB, IHC-P, ELISA with cross-reactivity to mouse, rat. Ideal for mitochondrial fatty acid oxidation research. Immunogen: recombinant fusion protein covering amino acids 31-210.

[Learn More](#)

Parameter	Value
Product Name	ACADL Rabbit pAb
Remarks/Alias	LCAD; ACAD4; ACADL
Species	Human
GeneID (Human)	33
GeneID	33
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 31-210 of human ACADL (NP_001599.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P28330
Protein Weight	48kDa
Shipping	Ice bag