

Anti-Acadm (Mcad) Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P11310

Item Number: CM0018780



Introduction

Affinity-purified rabbit polyclonal antibody against human ACADM (MCADH), raised against recombinant ACADM protein. Validated for WB, IF/ICC, ELISA; cross-reacts with mouse and rat. Ideal for studying fatty acid oxidation and mitochondrial disorders.

[Learn More](#)

Parameter	Value
Product Name	ACADM Rabbit pAb
Remarks/Alias	MCAD, ACAD1, MCADH
Species	Human
GeneID (Human)	34
GeneID	34
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human ACADM (NP_000007.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11310
Protein Weight	47kDa
Shipping	Ice bag