

Anti-Acadm Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P11310

Item Number: CM0009980



Introduction

Rabbit monoclonal antibody targeting human ACADM (MCAD) protein for mitochondrial fatty acid beta-oxidation research. Validated for WB, IHC-P, IF/ICC, ELISA; cross-reacts with human, mouse, rat. 47kDa.

[Learn More](#)

Parameter	Value
Product Name	ACADM Rabbit mAb
Remarks/Alias	MCAD; ACAD1; MCADH; ACADM
Species	Human
Gene ID (Human)	34
Gene ID	34
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11310
Protein Weight	47kDa
Shipping	Ice bag