


Anti-Acads Rabbit Polyclonal Antibody For Wb, Ip, Elisa - P16219

Item Number: CM0026764



Introduction

Rabbit polyclonal antibody against human ACADS (SCAD), cross-reactive with mouse and rat. Validated for WB, IP, and ELISA. Targets short-chain acyl-CoA dehydrogenase for mitochondrial fatty acid beta-oxidation research.

[Learn More](#)

Parameter	Value
Product Name	ACADS Rabbit pAb
Remarks/Alias	SCAD; ACAD3; ACADS
Species	Human
Gene ID (Human)	35
Gene ID	35
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human ACADS (NP_000008.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P16219
Protein Weight	44kDa
Shipping	Ice bag