


Anti-Acads / Scad Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P16219

Item Number: CM0003257



Introduction

Rabbit monoclonal antibody targeting human short-chain specific acyl-CoA dehydrogenase (SCAD). Validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Ideal for mitochondrial fatty acid beta-oxidation studies.

[Learn More](#)

Parameter	Value
Product Name	ACADS / SCAD Rabbit mAb
Remarks/Alias	SCAD; ACAD3; ACADS / SCAD
Species	Human
GeneID (human)	35
GeneID	35
Immunogen	Synthetic peptide (aa 313-412 of human ACADS/SCAD, P16219)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P16219
Protein Weight	44kDa
Shipping	Ice bag