

Anti-Acadvl Polyclonal Antibody For Wb, If/Icc, Elisa - P49748

Item Number: CM0029947



Introduction

Rabbit polyclonal antibody against human ACADVL, validated for WB, IF/ICC, ELISA. Detects endogenous ACADVL in human, mouse, and rat samples. Suitable for studying mitochondrial fatty acid beta-oxidation disorders.

[Learn More](#)

Parameter	Value
Product Name	ACADVL Rabbit pAb
Remarks/Alias	ACAD6; LCACD; VLCAD; ACADVL
Species	Human
GeneID (Human)	37
GeneID	37
Immunogen	Recombinant protein / Recombinant fusion protein containing a sequence corresponding to amino acids 1-140 of human ACADVL (NP_000009.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49748
Protein Weight	70kDa
Shipping	Ice bag