

Anti-Acap2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q15057

Item Number: CM0015169



Introduction

Rabbit polyclonal antibody against human ACAP2 (Arf-GAP with coiled-coil, ANK repeat and PH domain-containing protein 2), validated for WB, IF/ICC, and ELISA. Reacts with human and mouse. Predicted molecular weight 88 kDa.

[Learn More](#)

Parameter	Value
Product Name	ACAP2 Rabbit pAb
Remarks/Alias	CENTB2; CNT-B2; ACAP2
Species	Human
Gene ID (Human)	23527
Gene ID	23527
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 50-160 of human ACAP2 (NP_036419.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q15057
Protein Weight	88kDa
Shipping	Ice bag