


Anti-Acat1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P24752

Item Number: CM0003587



Introduction

Rabbit monoclonal antibody against human ACAT1, validated for WB, IHC-P, and ELISA. Reacts with human and rat samples. Ideal for beta-oxidation and ketone body metabolism research.

[Learn More](#)

Parameter	Value
Product Name	ACAT1 Rabbit mAb
Remarks/Alias	T2; MAT; ACAT; THIL; ACAT1
Species	Human
Gene ID (Human)	38
Gene ID	38
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 328-427 of human ACAT1 (P24752)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	P24752
Protein Weight	45kDa
Shipping	Ice bag