


Anti-Acat1 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P24752

Item Number: CM0007642



Introduction

High-quality rabbit monoclonal antibody targeting human ACAT1, a key mitochondrial enzyme in fatty acid beta-oxidation and ketone body metabolism. Validated for WB, IHC-P, IF/ICC, and ELISA; cross-reacts with human, mouse, and rat. Ideal for metabolic and mitochondrial research.

[Learn More](#)

Parameter	Value
Product Name	ACAT1 Rabbit mAb
Remarks/Alias	T2; MAT; ACAT; THIL; ACAT1
Species	Human
Gene ID (Human)	38
Gene ID	38
Immunogen	Recombinant protein corresponding to a sequence within amino acids 34-427 of human ACAT1 (NM_000019.4)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P24752
Protein Weight	45kDa
Shipping	Ice bag