


Anti-Acat1 Rabbit Polyclonal Antibody For Wb And Elisa - P24752

Item Number: CM0030511



Introduction

The Anti-ACAT1 Rabbit Polyclonal Antibody (P24752) is validated for western blot and ELISA, with reactivity across human, mouse, and rat. It targets mitochondrial acetyl-CoA acetyltransferase, a key enzyme in fatty acid beta-oxidation and ketone body metabolism.

[Learn More](#)

Parameter	Value
Product Name	ACAT1 Rabbit pAb
Remarks/Alias	T2, MAT, ACAT, THIL, ACAT1
Species	Human
Gene ID (Human)	38
Gene ID	38
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-145 of human ACAT1 (NP_000010.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P24752
Protein Weight	45kDa
Shipping	Ice bag