


Anti-Acat1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P24752

Item Number: CM0014984



Introduction

Rabbit polyclonal antibody against human ACAT1, suitable for WB, IHC-P, and ELISA. Reacts with human, mouse, and rat. Validated for research on fatty acid oxidation and ketone body metabolism.

[Learn More](#)

Parameter	Value
Product Name	ACAT1 Rabbit pAb
Remarks/Alias	T2, MAT, ACAT, THIL, ACAT1
Species	Human
GeneID (human)	38
GeneID	38
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-145 of human ACAT1 (NP_000010.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P24752
Protein Weight	45 kDa
Shipping	Ice bag