


Anti-Acat2 Rabbit Monoclonal Antibody For Wb, Elisa - Q9Bwd1

Item Number: CM0003751



Introduction

High-specificity rabbit monoclonal antibody targeting human ACAT2 (acetyl-CoA acetyltransferase 2). Validated for Western blot and ELISA; cross-reacts with mouse and rat. Ideal for cholesterol biosynthesis research.

[Learn More](#)

Parameter	Value
Product Name	ACAT2 Rabbit mAb
Remarks/Alias	ACAT2; acetyl-CoA acetyltransferase 2
Species	Human
Gene ID (Human)	39
Gene ID	39
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 150-250 of human ACAT2 (Q9BWD1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BWD1
Protein Weight	41kDa
Shipping	Ice bag