


Anti-Acat2 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Bwd1

Item Number: CM0024945



Introduction

Rabbit polyclonal antibody targeting human ACAT2, suitable for WB, IHC-P, IF/ICC, and ELISA. Detects mouse and rat. Cytosolic acetyl-CoA acetyltransferase involved in cholesterol biosynthesis. For IVD raw-material supply and research.

[Learn More](#)

Parameter	Value
Product Name	ACAT2 Rabbit pAb
Remarks/Alias	ACAT2
Species	Human
GeneID (Human)	39
Immunogen	Synthetic peptide within amino acids 159-277 of human ACAT2 (NP_005882.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BWD1
Protein Weight	41kDa
Shipping	Ice bag