


Anti-Acbd3 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9H3P7

Item Number: CM0021339



Introduction

High-affinity rabbit polyclonal antibody against human ACBD3 (GCP60) validated for WB, IF/ICC, and ELISA. Reacts with human, mouse, and rat. Ideal for Golgi apparatus and steroid biosynthesis research.

[Learn More](#)

Parameter	Value
Product Name	ACBD3 Rabbit pAb
Remarks/Alias	PAP7, GCP60, GOCAP1, GOLPH1, ACBD3
Species	Human
Gene ID (Human)	64746
Gene ID	64746
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-160 of human ACBD3 (NP_073572.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H3P7
Protein Weight	61kDa
Shipping	Ice bag