


Anti-Acc1 Monoclonal Antibody For Wb, Ihc-P, Ip, Elisa - Q13085

Item Number: CM0002149



Introduction

Rabbit monoclonal antibody targeting human ACC1, a key enzyme in fatty acid biosynthesis. Suitable for WB, IHC-P, IP, ELISA, with cross-reactivity to mouse and rat. Ideal for metabolic research.

[Learn More](#)

Parameter	Value
Product Name	ACC1 Rabbit mAb
Remarks/Alias	ACC; ACAC; ACC1; ACCA; Acac1; hACC1; ACACAD; ACCalpha; ACACalpha
Species	Human
GeneID (human)	31
GeneID	31
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 350-450 of human ACC1 (Q13085).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13085
Protein Weight	266kDa
Shipping	Ice bag