


## Anti-Acc1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - Q13085

Item Number: CM0015586



### Introduction

Rabbit polyclonal antibody targeting human ACC1 (Acetyl-CoA carboxylase 1), suitable for WB, IHC-P, IF/ICC, IP, ELISA. Detects human, mouse, rat. For fatty acid biosynthesis studies.

[Learn More](#)

| Parameter        | Value                                                                                                               |
|------------------|---------------------------------------------------------------------------------------------------------------------|
| Product Name     | ACC1 Rabbit pAb                                                                                                     |
| Remarks/Alias    | ACC; ACAC; ACC1; ACCA; Acac1; hACC1; ACACAD; ACCalpha; ACACalpha                                                    |
| Species          | Human                                                                                                               |
| GeneID (Human)   | 31                                                                                                                  |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 2266-2346 of human ACC1 (NP_942133.1) |
| Source           | Rabbit                                                                                                              |
| Category         | Polyclonal Antibodies                                                                                               |
| Application      | WB, IHC-P, IF/ICC, IP, ELISA                                                                                        |
| Cross Reactivity | Human, Mouse, Rat                                                                                                   |
| SWISS            | Q13085                                                                                                              |
| Protein Weight   | 266 kDa                                                                                                             |
| Shipping         | Ice bag                                                                                                             |