


Anti-Acd Rabbit Polyclonal Antibody For Wb And Elisa - Q96Ap0

Item Number: CM0017125



Introduction

Anti-ACD Rabbit Polyclonal Antibody, validated for WB and ELISA with human samples. Targets the shelterin protein ACD (PIP1/PTOP/TPP1) involved in telomere maintenance and protection.

[Learn More](#)

Parameter	Value
Product Name	ACD Rabbit pAb
Remarks/Alias	PIP1; PTOPI; TPP1; TINT1; ACD
Species	Human
GeneID (human)	65057
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human ACD (NP_001075955.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q96AP0
Protein Weight	49kDa
Shipping	Ice bag