


Anti-Ace1 Polyclonal Antibody - P12821

Item Number: CM0028226



Introduction

High-quality rabbit polyclonal antibody against human ACE1 (angiotensin-converting enzyme) for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for cardiovascular and neurobiology research. Protein weight: 150kDa.

[Learn More](#)

Parameter	Value
Product Name	ACE1 Rabbit pAb
Remarks/Alias	DCP; ACE1; DCP1; CD143
Species	Human
Gene ID (Human)	1636
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 554-758 of human ACE1 (NP_000780.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P12821
Protein Weight	150kDa
Shipping	Ice bag