


Anti-Ace1 Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P12821

Item Number: CM0005417



Introduction

High-quality Anti-ACE1 rabbit monoclonal antibody validated for WB, IF-P, IHC-P, ELISA. Detects human, mouse, rat ACE1 (P12821). Ideal for cardiovascular, neuropeptide, and hematopoietic stem cell research.

[Learn More](#)

Parameter	Specification
Product Name	ACE1 Rabbit mAb
Remarks/Alias	DCP; ACE1; DCP1; CD143
Species	Human
GeneID (Human)	1636
GeneID	1636
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-1261 of human ACE1 (NP_000780.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P12821
Protein Weight	150 kDa
Shipping	Ice bag