

Anti-Ace2 Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q9Byf1

Item Number: CM0010087



Introduction

Rabbit monoclonal ACE2 antibody validated for WB, IHC-P, IF-P, and ELISA. Reacts with human, mouse, and rat. Targets the 92 kDa angiotensin-converting enzyme 2, a key regulator of cardiovascular homeostasis and SARS-CoV-2 receptor. Ideal for renin-angiotensin system and coronavirus research.

[Learn More](#)

Parameter	Value
Product Name	ACE2 Rabbit mAb
Remarks/Alias	ACEH; ACE2
Species	Human
GeneID (Human)	59272
GeneID	59272
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-740 of human ACE2 (NP_068576.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BYF1
Protein Weight	92 kDa
Shipping	Ice bag