

Anti-Ace2 Rabbit Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q9Byf1

Item Number: CM0016961



Introduction

Polyclonal antibody targeting human ACE2, validated for Western blot, immunofluorescence, immunohistochemistry, and ELISA. Recognizes human, mouse, and rat ACE2. Suitable for renin-angiotensin system and SARS-CoV-2 receptor studies.

[Learn More](#)

Parameter	Specification
Product Name	ACE2 Rabbit pAb
Remarks/Alias	ACEH; ACE2
Species	Human
GeneID (Human)	59272
GeneID	59272
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 18-740 of human ACE2 (NP_068576.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BYF1
Protein Weight	92 kDa
Shipping	Ice bag