

## Anti-Ace2 Rabbit Polymab® For Wb, If-P, Ihc-P, Elisa - Q9Byf1

Item Number: CM0002738



### Introduction

CamelBio supplies Anti-ACE2 Rabbit PolymAb® (Q9BYF1), validated in WB, IF-P, IHC-P, and ELISA. Reacts with human, mouse, and rat ACE2, a key enzyme of the renin-angiotensin system and receptor for SARS coronaviruses. Ideal for cardiovascular and infectious disease research.

[Learn More](#)

| Parameter        | Details                |
|------------------|------------------------|
| Product Name     | ACE2 Rabbit PolymAb®   |
| Remarks/Alias    | ACEH                   |
| Species          | Human                  |
| GeneID           | 59272                  |
| Immunogen        | Synthetic peptide      |
| Source           | Rabbit                 |
| Category         | Monoclonal Antibodies  |
| Application      | WB, IF-P, IHC-P, ELISA |
| Cross Reactivity | Human, Mouse, Rat      |
| SWISS            | Q9BYF1                 |
| Protein Weight   | 53 kDa/92 kDa          |
| Shipping         | Ice bag                |