

## Anti-Acha4 Rabbit Monoclonal Antibody For Wb, Elisa - P43681

Item Number: CM0003872



### Introduction

Rabbit monoclonal antibody against human ACHA4 (CHRNA4), validated for Western blot and ELISA. Detects endogenous levels of ACHA4 in human, mouse, and rat samples. Ideal for neuroscience research.

[Learn More](#)

| Parameter        | Value                                                                                               |
|------------------|-----------------------------------------------------------------------------------------------------|
| Product Name     | ACHA4 Rabbit mAb                                                                                    |
| Remarks/Alias    | EBN; BFNC; EBN1; NACHR; NACRA4; NACHRA4; ACHA4                                                      |
| Species          | Human                                                                                               |
| GeneID (Human)   | 1137                                                                                                |
| GeneID           | 1137                                                                                                |
| Immunogen        | A synthetic peptide corresponding to a sequence within amino acids 500-600 of human ACHA4 (P43681). |
| Source           | Rabbit                                                                                              |
| Category         | Monoclonal Antibodies                                                                               |
| Application      | WB, ELISA                                                                                           |
| Cross Reactivity | Human, Mouse, Rat                                                                                   |
| SWISS            | P43681                                                                                              |
| Protein Weight   | 70kDa                                                                                               |
| Shipping         | Ice bag                                                                                             |