


Anti-Ache Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P22303

Item Number: CM0029485



Introduction

Rabbit polyclonal antibody against human acetylcholinesterase (ACHE, P22303). Validated for WB, IF/ICC, ELISA; cross-reacts with mouse and rat. Targets the C-terminal region for detecting membrane-bound and soluble forms.

[Learn More](#)

Parameter	Value
Product Name	ACHE Rabbit pAb
Remarks/Alias	YT; ACEE; ARACHE; N-ACHE; ACHE
Species	Human
GeneID (human)	43
GeneID	43
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 398-617 of human ACHE (NP_056646.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P22303
Protein Weight	68kDa
Shipping	Ice bag