


Anti-Ackr3 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P25106

Item Number: CM0018301



Introduction

Anti-ACKR3 Rabbit Polyclonal Antibody for WB, IF/ICC, ELISA. Detects human, mouse, rat ACKR3 (CXCR7/RDC1), an atypical chemokine receptor regulating chemokine levels and MAPK pathway.

[Learn More](#)

Parameter	Value
Product Name	ACKR3 Rabbit pAb
Remarks/Alias	RDC1; CXCR7; RDC-1; CMKOR1; CXC-R7; CXCR-7; GPR159; ACKR3
Species	Human
Gene ID (Human)	57007
Gene ID	57007
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 150-250 of human ACKR3 (NP_064707.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P25106
Protein Weight	41kDa
Shipping	Ice bag