


Anti-Acly Polyclonal Antibody For Wb, Ihc-P, If, Elisa - P53396

Item Number: CM0017545



Introduction

Rabbit polyclonal antibody against human ACLY (ATP-citrate synthase). Validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Ideal for metabolism research and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	ACLY Rabbit pAb
Remarks/Alias	ACL; ATPCL; CLATP; ACLY
Species	Human
GeneID (Human)	47
GeneID	47
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 752-1101 of human ACLY (NP_001087.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P53396
Protein Weight	121kDa
Shipping	Ice bag