


Anti-Acly Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P53396

Item Number: CM0005047



Introduction

Anti-ACLY rabbit monoclonal antibody validated for Western blot, IHC-P, and ELISA. Reacts with human, mouse, and rat ACLY. Suitable for lipid metabolism, cancer metabolism, and metabolic disease research.

[Learn More](#)

Parameter	Description
Product Name	ACLY Rabbit mAb
Remarks/Alias	ACL; ATPCL; CLATP; ACLY
Species	Human
GeneID (Human)	47
GeneID	47
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 752-1101 of human ACLY (NP_001087.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P53396
Protein Weight	121kDa
Shipping	Ice bag