

Anti-Acmsd Rabbit Polyclonal Antibody For Wb, Elisa - Q8Tdx5

Item Number: CM0023578



Introduction

Rabbit polyclonal antibody against human ACMSD, suitable for Western blot and ELISA. Cross-reacts with mouse and rat. Validated for tryptophan metabolism research. Predicted molecular weight 38 kDa.

[Learn More](#)

Parameter	Value
Product Name	ACMSD Rabbit pAb
Remarks/Alias	ACMSD
Species	Human
GeneID (Human)	130013
GeneID	130013
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-336 of human ACMSD (NP_612199.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q8TDX5
Protein Weight	38kDa
Shipping	Ice bag