


## Anti-Aconitase 1 (Aco1) Polyclonal Antibody For Wb, Elisa - P21399

Item Number: CM0026701



### Introduction

Rabbit polyclonal antibody targeting human Aconitase 1 (ACO1/IRP1) for WB and ELISA. Validated in human, mouse, and rat samples. Essential for iron metabolism and TCA cycle research.

[Learn More](#)

| Parameter        | Value                                                                                                         |
|------------------|---------------------------------------------------------------------------------------------------------------|
| Product Name     | Aconitase 1 (ACO1) Rabbit pAb                                                                                 |
| Remarks/Alias    | IRP1; ACONS; HEL60; IREB1; IREBP; IREBP1; Aconitase 1 (ACO1)                                                  |
| Species          | Human                                                                                                         |
| GeneID (human)   | 48                                                                                                            |
| GeneID           | 48                                                                                                            |
| Immunogen        | Synthetic peptide corresponding to a sequence within amino acids 100-200 of human Aconitase 1 (ACO1) (P21399) |
| Source           | Rabbit                                                                                                        |
| Category         | Polyclonal Antibodies                                                                                         |
| Application      | WB, ELISA                                                                                                     |
| Cross Reactivity | Mouse, Rat                                                                                                    |
| SWISS            | P21399                                                                                                        |
| Protein Weight   | 98kDa                                                                                                         |
| Shipping         | Ice bag                                                                                                       |