


Anti-Aconitase 1 (Aco1) Polyclonal Antibody For Wb - P21399

Item Number: CM0025509



Introduction

Rabbit polyclonal antibody against human Aconitase 1 (ACO1, also IRP1) for Western blot and ELISA. Validated in human, mouse, and rat. Detects a 98kDa protein.

[Learn More](#)

Parameter	Value
Product Name	Aconitase 1 (ACO1) Rabbit pAb
Remarks/Alias	IRP1; ACONS; HEL60; IREB1; IREBP; IREBP1; Aconitase 1 (ACO1)
Species	Human
GeneID (Human)	48
GeneID	48
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-170 of human Aconitase 1 (ACO1) (NP_002188.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P21399
Protein Weight	98kDa
Shipping	Ice bag