


Anti-Aconitase 1 (Aco1) Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P21399

Item Number: CM0011250



Introduction

Rabbit monoclonal antibody against human Aconitase 1 (ACO1/IRP1), validated for WB, IF/ICC, and ELISA. Reacts with human, mouse, and rat. Suitable for iron metabolism and gene regulation research.

[Learn More](#)

Parameter	Value
Product Name	Aconitase 1 (ACO1) Rabbit mAb
Remarks/Alias	IRP1; ACONS; HEL60; IREB1; IREBP; IREBP1; Aconitase 1 (ACO1)
Species	Human
GeneID (Human)	48
GeneID	48
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human Aconitase 1 (ACO1) (P21399)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P21399
Protein Weight	98kDa
Shipping	Ice bag