

Anti-Aconitase 2 (Aco2) Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q99798

Item Number: CM0006146

Introduction

Rabbit monoclonal antibody targeting human Aconitase 2 (ACO2), validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for studying mitochondrial TCA cycle enzyme function and related neurodegenerative disorders.

[Learn More](#)

Parameter	Value
Product Name	Aconitase 2 (ACO2) Rabbit mAb
Remarks/Alias	ICRD; OCA8; OPA9; ACONM; HEL-S-284; Aconitase 2 (ACO2)
Species	Human
GeneID (Human)	50
GeneID	50
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of human Aconitase 2 (ACO2) (Q99798).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99798
Protein Weight	85kDa
Shipping	Ice bag