

Anti-Aconitase 2 (Aco2) Rabbit Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - Q99798

Item Number: CM0026798

Introduction

Rabbit polyclonal antibody against human Aconitase 2 (ACO2), validated for WB, IF/ICC, IHC-P, ELISA. Cross-reacts with mouse and rat. Suitable for mitochondrial metabolism and TCA cycle research.

[Learn More](#)

Parameter	Value
Product Name	Aconitase 2 (ACO2) Rabbit pAb
Remarks/Alias	ICRD; OCA8; OPA9; ACONM; HEL-S-284; Aconitase 2 (ACO2)
Species	Human
Gene ID (Human)	50
Gene ID	50
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 501-780 of human Aconitase 2 (Aconitase 2 (ACO2)) (NP_001089.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99798
Protein Weight	85kDa
Shipping	Ice bag